



sonosim®

SonoSim LiveScan®
User Guide

Table of Contents


Using SonoSim LiveScan - Getting Started	3
Launch SonoSim LiveScan	3
Tag Setup	4
Calibrate the Probe	5
Choose Your Tags	6
Apply Your Tags	7
Save Setup	10
Edit Tag Setup	11
Select a Case	12
Start Scanning	13
SonoSim LiveScan Features	14
Cleaning & Disinfection	15
Member Support & Personal Onboarding	16
Regulatory Compliance	17

Using SonoSim LiveScan - Getting Started

To begin using SonoSim LiveScan, please follow these steps:

1- Launch SonoSim LiveScan

From the Main Menu, select the button "Apply" that displays the SonoSim LiveScan logo.




2- Tag Setup

SonoSim LiveScan tags must be precisely positioned to ensure accurate ultrasound window representation. Follow Tag Setup instructions on screen for each tag you wish to place. From the SonoSim LiveScan Tag Setup Menu you will have three options to choose from:

New: If you have not already set up tags for a new manikin or human volunteer select the 'New' button. You will be walked through the tag setup for each tag. Once completed, you can save your tag setup for that specific manikin or human volunteer.

Saved: Load saved setup configurations you have previously created by clicking the 'Saved' button. Once you have pressed the "Saved" button you will be brought to the list of all your saved setups.

Default (Not Recommended): If you wish to bypass tag setup, you can use the default configurations. *Due to the unique shape of each model, ultrasound images in the Default Setup may not match your probe positioning. Creating a new tag setup ensures that the ultrasound images display correctly on your model and gives an optimal experience.*




3- Calibrate the SonoSim LiveScan Probe

Step 1: Hold the Probe upright on a flat surface along the sagittal plane of the patient.

Step 2: Point the indicator cephalad.

Step 3: Press "Calibrate" or "C" key.

SonoSim® LiveScan Probe Calibration



Place the Probe in **PATIENT POSITION**:

1. Hold the Probe upright on a flat surface along the sagittal plane of a patient
2. Point the indicator cephalad


Calibrate

Note: SonoSim LiveScan tags must NOT be used on anyone with an electronic pacemaker, implantable cardioverter-defibrillator, or other implantable electronic medical device.


Using SonoSim LiveScan

4- Choose your SonoSim LiveScan Tags

Choose either semi-permanent or disposable SonoSim LiveScan tags, depending on whether the tags will be applied to a manikin or human volunteer. The disposable tags (silver) have a hypoallergenic adhesive and are specifically designed for single use application onto human skin. The semi-permanent tags (blue) are designed for long-term application onto manikin.



Semi-permanent manikin SonoSim LiveScan tags



Disposable human SonoSim LiveScan tags

Note: SonoSim LiveScan tags must NOT be used on anyone with an electronic pacemaker, implantable cardioverter-defibrillator, or other implantable electronic medical device.

5- Apply your SonoSim LiveScan Tags

Step 1: Begin by placing the first tag on your manikin or volunteer, according to the instructions on screen.

Step 2: Place the probe flat on the SonoSim LiveScan Tag and click 'Align'.

Tag Setup

Standard Tags | OB/GYN Tags


Left Chest	Apical	RUQ	Proximal IVC	X-7	Left Groin	Suprapubic	Left Eye
Right Chest	Parasternal	Subcostal	LUQ	Mid Aorta	Right Groin	Head	Right Eye

Step 1. Place Left Chest tag on model and scan.
Apply tag along the left midclavicular line between the second to fourth rib interspaces.
Apply superior to nipple line.

Step 2. Place probe flat on Left Chest tag and click Align.

Align

Restart Setup | Finish Setup



Using SonoSim LiveScan

Apply and setup each tag you plan to use. It is not required to set up all tags available.

After you complete each tag setup, the image of the tag on the screen will fill in blue to indicate you have aligned the tag to your manikin or volunteer.

Tag Setup

Standard Tags | OB/GYN Tags

Left Chest	Apical	RUQ	Proximal IVC	X-7	Left Groin	Suprapubic	Left Eye
Right Chest	Parasternal	Subcostal	LUQ	Mid Aorta	Right Groin	Head	Right Eye

Step 1. Place Mid Aorta tag on model and scan.
Apply tag three fingerbreadths inferior to the Proximal IVC tag. Tag should be approximately midway between umbilicus and xiphoid process in an adult model.

Step 2. Place probe flat on Mid Aorta tag and click Align.

Align

Restart Setup | Finish Setup

Using SonoSim LiveScan

If you need to place OB/GYN tags, select the 'OB/GYN Tags' tab.

Note: After setting up Standard Tags you will automatically be directed to OB/GYN tag setup. It is not required to setup all tags - only the ones you wish to use. Skip any unwanted tags by selected the next tag you wish to set up.

Tag Setup

Standard Tags **OB/GYN Tags**

Lower Uterus Infraumbilical Para Ut RLQ **Para Ut LLQ** Para Ut LMQ



Upper Uterus Supraumbilical Para Ut RUQ Para Ut LUQ Para Ut RMQ

Step 1. Place Para Ut LLQ tag on model and scan.
Apply tag in the left lower quadrant five fingerbreadths inferior lateral from the umbilicus. Tag should be below the umbilicus.

Step 2. Place probe flat on Para Ut LLQ tag and click Align.



Align

Restart Setup Finish Setup




6- Save Setup


Once you have placed and setup all of the tags that you will use for your simulation, click "Finish Setup" and the simulator will prompt you to name and save your template.



Once you have entered the SonoSim LiveScan Simulator, you will see the name of the setup in use, in the lower left corner of the screen.



7- Edit Tag Setup

If you wish to create a new setup (A) or edit your current tag setup (B), navigate to the  button (C) on the bottom right side of the simulator.




From the SonoSim LiveScan Tag Setup menu you can also load and edit configurations you have previously created by clicking the “Saved” button. Once you have pressed the “Saved” button you will be brought to the list of all your saved setups.

Once you are in the saved set up menu, hovering over each setup will give you the option to “Load”, “Edit”, or “Delete” any setup.

“Load” will bring you directly to the simulator with the setup you selected.

“Edit” will allow you to re-align the tags in the selected setup.



“Delete” will allow you to delete the selected setup.



Using SonoSim LiveScan

8- Select a case

Select the  button and select a SonoSim LiveScan Category.



9- Start scanning

Place the SonoSim LiveScan Probe on a tag and start scanning. To change the image window, simply lift and position the probe on another SonoSim LiveScan tag.

If your images do not show as expected when you put the probe down, modify your tag setup ([see Edit Tag Setup instructions](#)).



SonoSim LiveScan Features

The image shows a screenshot of the SonoSim LiveScan software interface. The central area displays a B-mode ultrasound image. Surrounding the image are several control panels and toolbars. Callout boxes with blue borders and white text point to specific features:

- Perform various measurements:** Points to the 'Caliper' and 'Worksheet' icons in the top-left toolbar.
- Adjust image depth:** Points to the 'Depth' control in the top-right vertical toolbar.
- Adjust image brightness:** Points to the 'Gain' control in the top-right vertical toolbar.
- Display available Doppler:** Points to the 'Doppler' icon in the top-right vertical toolbar.
- Flip the Probe indicator:** Points to the 'Orientation' icon in the top-right vertical toolbar.
- Freeze, save, review, & annotate images & data:** Points to the 'STORAGE' section containing 'Freeze', 'Save', 'Review', and 'Annotate' icons in the bottom-right vertical toolbar.
- Navigate to available imaging view:** Points to a row of six circular icons at the bottom center, with the fourth icon highlighted in blue.
- Exit to the Main Menu:** Points to the back arrow icon at the bottom of the right-hand vertical toolbar.
- Navigate to next case:** Points to the left and right arrow icons in the right-hand vertical toolbar.
- Case history:** Points to the 'Hx' icon in the right-hand vertical toolbar.
- Lock Probe:** Points to the lock icon in the right-hand vertical toolbar.
- Image interpretation by an expert instructor:** Points to the play button icon in the right-hand vertical toolbar.
- Compress a target vessel:** Points to the vessel compression icon in the right-hand vertical toolbar.
- Calibrate Probe for optimal scanning experience:** Points to the crosshair icon in the right-hand vertical toolbar.
- Return to the tag setup menu, edit your tags, or access the user guide:** Points to the question mark icon in the right-hand vertical toolbar.

Cleaning & Disinfection

To clean your SonoSim LiveScan Probe:

1. Damp a small piece of dry cloth with 70% Isopropyl Alcohol (IPA)
2. Wipe the exterior of the Probe gently, ensuring no liquid enters the seams

Caution: Never submerge the probe in liquid.

To clean your SonoSim LiveScan Tags:

1. Damp a small piece of dry cloth with 70% Isopropyl Alcohol (IPA)
2. Wipe the top of the tags gently

Caution: Never apply any cleaning agent to the adhesive side of the tags.



Member Support & Personal Onboarding

Online Resources

For online service and support information, please visit sonosim.com/support and watch tutorials at sonosim.com/videos

SonoSim Member Support

For a personal demonstration or technical help using the features, please contact membersuccess@sonosim.com or call us at 855.873.7666 (M-F 6am-6pm, Sat 8am-4pm PT). Enjoy your ultrasound training journey!

Federal Communications Commission Statement

This device complies with Part 15 of the FCC Rules and RSS-210 of IC rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna
- Increase the separation between the equipment and the receiver
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected
- Consult the dealer or an experienced radio/TV technician for help.

The term "IC" before the equipment certification number only signifies that the Industry Canada technical specifications were met. Changes or modification not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Déclaration de la Commission Fédérale des Communications

Cet appareil est conforme à la section 15 des règles de la FCC et RSS-210 des règles IC. Son fonctionnement est sujet aux deux conditions suivantes :

1. Cet appareil ne doit pas causer des interférences nuisibles, et
 2. Cet appareil doit accepter toute interférence reçue, y compris les interférences qui peuvent provoquer un fonctionnement indésirable.
- Cet équipement a été testé et jugé conforme aux limites d'un appareil numérique de classe B, conformément à la section 15 des règles de la FCC. Ces limites sont conçues pour fournir une protection raisonnable contre les interférences nuisibles dans une installation résidentielle. Cet équipement génère, utilise et peut émettre une énergie radiofréquence et, s'il n'est pas installé et utilisé conformément aux instructions, peut causer des interférences nuisibles aux communications radio. Cependant, il n'y a aucune garantie que des interférences ne se produiront pas dans une installation particulière. Si cet équipement provoque des interférences nuisibles à la réception radio ou télévision, ce qui peut être déterminé en allumant et en éteignant l'équipement, l'utilisateur est encouragé à essayer de corriger l'interférence par une ou plusieurs des mesures suivantes :
- Réorienter ou déplacer l'antenne réceptrice.
 - Augmenter la distance entre l'équipement et le récepteur.
 - Connecter l'équipement à une prise sur un circuit différent de celui sur lequel le récepteur est branché.
 - Consulter le revendeur ou un technicien radio/TV expérimenté pour obtenir de l'aide.

Le terme "IC" avant le numéro de certification de l'équipement signifie uniquement que les spécifications techniques d'Industrie Canada ont été respectées.

Les changements ou modifications non expressément approuvés par la partie responsable de conformité pourraient annuler l'autorité de l'utilisateur à utiliser l'équipement .

Waste Handling



Dispose of in accordance with local requirements and regulations.

Certification



The product is in compliance with the essential requirements of Council Directive 1999/5/EC on radio and telecommunications terminal equipment (R&TTE). The product is in compliance with Council Directive 2011/65/EU on restriction of the use of certain hazardous substances (RoHS).